


## Procedure/Surgical Masks 3-Ply Adult Level 3



### Description & Classification

Model	3-Ply Adult, Level 3
Product Name	Procedure/Surgical Mask, Non-Sterile
Strap Type	Earloop
Size	Standard, Fits most adults
Sample	
Measurements	W: 6.9" x H: 3.7"
Colour	Blue, Black, White
Manufactured In	Canada
Use Range	Medical Facilities Personal Protection

### Material and Performance Specification

Material	Non-woven Fabric
Contains Latex	No
Bacterial Filter Efficiency	≥99.88%
Differential Pressure	≤3.17
Particulate Filter Efficiency	≥99.7%
Fluid Resistance	≥160mm Hg
Flammability Rating	Class I

### Performance Requirements

Test	ASTM F2100			3-Ply Adult, Level 3
	Level 1	Level 2	Level 3	
Bacterial Filtration Efficiency	≥95%	≥98%	≥98%	Pass
Differential Pressure	< 5 mmH <sub>2</sub> O	< 6 mmH <sub>2</sub> O	< 6 mmH <sub>2</sub> O	Pass
Particulate Filtration Efficiency	≥95%	≥98%	≥98%	Pass
Fluid Resistance	80 mm Hg	120 mm Hg	160 mm Hg	Pass
Flammability Rating	Class 1	Class 1	Class 1	Pass



### Packaging Specification

Mask Weight	0.11 oz
Masks in Box	50
Weight of Box with Masks	0.41 lbs
Box Dimensions	L: 7.5" x W: 4" x H: 4"
Masks Per Case	1000
Boxes Per Case	20
Weight of Case	9.5 lbs
Case Dimensions	L: 20" x W: 15.5" x H: 8.75"
Masks Per Pallet	42000
Cases Per Pallet	42
Pallet Dimensions	L: 41" x W: 48" x H: 66"