

# SuperVista

## USB IntraOral Camera System



- Superior Image Quality
- Take dental diagnostics to the next level with this easy –to-use imaging tool
- The convenience of an intraoral camera that integrates with a digital dental radiography system
- Seamlessly compatible with the **EVASoft** Digital Dental Imaging System, the imaging software used by EVA digital X-ray dental offices
- USB 2.0 Connectivity—No external power supplies or proprietary capture cards needed



[www.ez2000digital.com](http://www.ez2000digital.com)

EZ 2000, Inc.  
Digital Dental Imaging Specialists  
800-273-5033

## EZ 2000 Dental Imaging Solutions for Every Dental Practice

The **SuperVista** intraoral camera is an easy-to-use hand-held USB camera. This wand-style camera uses digital photography for diagnostic evaluation and client or patient presentation. Simply position the camera close to a region of interest and use the single-click shutter button to capture as many images as you would like. After a picture is captured, it saves in the software as the camera returns to "live feed" video mode, ready for the next picture.

Many dentists use **EVA** digital X-ray sensors with the **EVASoft** imaging software. Using **EVASoft** with the **Supervista** camera means that you can save dental photographs alongside dental radiographs in a patient's file. This expansion of **EVASoft** usage provides more information to referring dentists and patients, thus allowing an office quick explanations through visual records.

The **SuperVista** is widely regarded as a critical aid in documentation and explanation.

### Each SuperVista Camera System Includes:

- 1 USB Camera
- 1 USB Cable
- 1 Wall/Table mounting holster
- 2 Mounting screws
- 1 Driver CD
- 1 Pkg. of protective sheaths
- 1 User Manual



Camera



Holster

### Specifications:

Working distance:	5~35mm
Light Output:	>=500 Lux
Resolution:	640 X 480 ~ 1280 X 1024
Field of View:	90 deg. +- 5%
Power:	5Vdc (via USB 2.0)
Video output:	USB 2.0 (full motion digital)
Imager:	1/4 inch CMOS
Illumination:	(6) LEDs
Temp & Humidity:	5 deg. ~ 68deg. C / <= 80%
Handpiece weight:	40 grams
Operating System:	Windows 7/Vista/XP
Drivers:	Windows 7: 64 bit & 32 bit Windows Vista & XP; 32 bit



EZ 2000, Inc.  
1801 Century Park East  
24th Floor  
Los Angeles, CA 90067

tel: 310-204-6134  
fax: 310-204-1440  
[www.ez2000digital.com](http://www.ez2000digital.com)